



Providing Excellence in Science Education

IPC Irwin Low Tension PSU (20V)

Part Number IPC-0195-P

User Instructions



INTRODUCTION.

The 20V Variable Low Tension (L.T.) PSU is a heavy-duty power supply designed to provide both continuously variable AC and unsmoothed, full-wave rectified, DC voltage supplies up to a current of 6 Amps, or 8 Amps for short periods. The high current capability makes the PSU ideal for operating motors, lamps, heating elements etc. plus demonstrating electromagnetism and performing magnetic experiments. Please note that the PSU now comes with a carrying handle fitted to the lid.

DESCRIPTION.

The Power Supply provides both continuously variable AC and DC voltages up to 20V, both outputs may be used simultaneously up to a combined output current of 6 Amps, or 8 Amps for short periods.

The control panel incorporates a voltage adjustment control knob, which is used to set the required output voltage, the approximate voltages are shown surrounding the control. The control panel also incorporates four 4mm output sockets, two White sockets for the AC output and a Red and Black socket for the DC output, a resettable overload cut-out switch is also provided on the panel.

The rear panel incorporates an illuminated (green) mains on/off switch and a mains input IEC socket containing an in-built protection fuse. Normally the unit is supplied with a separate UK plug to IEC plug, 2.0m, mains lead.

IPC Irwin

Holker School, Cark in Cartmel, Grange over Sands, Cumbria LA11 7PQ, United Kingdom

Call: 015395 58555 **Email:** sales@ipcirwin.com



Providing Excellence in Science Education

INSTRUCTIONS.

To set the required AC or DC output voltage simply rotate the voltage control knob to the desired voltage using the graduated voltage scale, then connect the power leads to the desired voltage output sockets, ie. the two white 4mm sockets for AC, and the red (positive) and black (negative) for DC.

If required both the AC and DC outputs may be used simultaneously up to a combined output current of 6 Amps, or 8 Amps for short periods. Please note, when using the PSU in this configuration you will only be able to set the AC and DC voltages to the same level.

The cut-out switch operates if there is a short circuit or overload on either of the outputs in excess of the maximum 8 Amps output. Should the reset button pop out, switch off the supply, rectify the fault, switch on again and press the button in. Please note that the button takes about 5 seconds to reset and will not work if pushed in within this time.

Should the mains on/off switch fail to illuminate and there is no output power from the unit, then first check the 2.5A fuse at the rear of the unit (located in the mains IEC socket), if this is found to be intact then check the 3A fuse in the mains plug (if UK plug).

Please note:- when there is no load connected (or a very small load) the AC and DC output voltages are higher than shown on the voltage control scale, for a more accurate measurement of the output voltage a suitable voltmeter or multimeter should be employed.

SPECIFICATIONS.

Electrical Supply	= 220-240VAC, 50-60Hz
Fuse Ratings	= IEC Socket Inlet: 2.5A, UK Mains Plug: 3A
Dimensions	= 255 x 230 x 120mm overall
Mass	= 4.65kg
AC Output Voltage Range	= 0 to 20VAC
AC Output Current	= 6A (8A for a short periods)
DC Output Voltage Range	= 0 to 20VDC
DC Output Current	= 6A (8A for a short periods)

Please Note: This Equipment is for Education, Training, or research Purposes Only.

IPC Irwin

Holker School, Cark in Cartmel, Grange over Sands, Cumbria LA11 7PQ, United Kingdom

Call: 015395 58555 **Email:** sales@ipcirwin.com



Providing Excellence in Science Education

IMPORTANT- ELECTRICAL SAFETY NOTICE

Only connect the product to single phase mains electrical supply (with neutral nominally at each potential), also check that the supply voltage and frequency are within the products stated range. Before connecting to the mains electrical supply, examine the information on the products rating label (the value of the fuse fitted to the product is also stated), ensure that the mains plug or outlet circuit is fitted with an appropriate fuse of higher value.

WARNING: THIS PRODUCT MUST BE EARTHED

This product is designed for Electrical Safety Class 1.

If a mains lead is not fitted or supplied, connect the wires to a Standard UK, 3-pin plug as follows:-

Green/Yellow wire to the terminal marked: E (Earth) or coloured Green or Green/Yellow.

Brown wire to terminal marked: L (Live) or coloured Brown.

Blue wire to terminal marked: N (Neutral) or coloured Blue.

LIVE PARTS SHOULD NEVER BE EXPOSED UNLESS THE PRODUCT HAS BEEN SWITCHED OFF AND ISOLATED FROM THE MAINS ELECTRICITY SUPPLY.

Thank you for choosing this IPC Irwin product, which has specifically been designed for education and we know it will bring you good service long into the future. This is just one product of a wide range of UK designed and manufactured products to enhance a STEM education, please see www.ipcirwin.com website for more details on the products and service we provide to support teaching and learning in science and technology.

IPC Irwin also provide Portable Appliance Testers and a Calibration service, for more details please go to www.ipcirwin.com or call our sales and technical support line on 015395 58555.



IPC Irwin

Holker School, Cark in Cartmel, Grange over Sands, Cumbria LA11 7PQ, United Kingdom

Call: 015395 58555 **Email:** sales@ipcirwin.com



Providing Excellence in Science Education

Declaration of Conformity (In Accordance with EN ISO 17050-1:2010).

Manufacturer: **IPC Irwin, Holker School, Cark-in-Cartmel, Grange-over-Sands, Cumbria, LA11 7PQ, United Kingdom.**

Product/Part No: **Low Tension (20V) PSU (230V AC, 50Hz main voltage)/IPC-0195-P**

This conformity certificate approves the compliance of the product with the essential safety requirements of the following UK Standards and Regulations:

SI 2005 No.1803	General Product Safety Regulation (implements Directive 2001/95/EC).
SI 2016/1101	Electrical Equipment (Safety) Regulation (implements Directive 2014 /35/EU.
EN 61010-1: 2010	Electrical Requirements for Laboratory Test & Measurement Equipment.
EN 61000-6-1: 2019	EMC, Generic Standards, Immunity for Residential, Commercial and Light Industrial Environments.
EN 61000-6-3: 2019	EMC, Generic Standards, Emission for Residential, Commercial and Light Industrial Environments.
SI 2012 No. 3032	Restriction on Hazardous Substances (RoHS) in Electrical and Electronic Equipment Regulation, (implements Directive 2011/65/EU).



Mark can be used only in the case of conformity assessment according to all relevant UK Standards and Regulations.

EU Declaration of Conformity (In accordance with European Parliament and Council Decision No: 768/2008/EC Annex III).

This conformity certificate approves the compliance of the product with the essential safety requirements of the following EC/EU Directives (Manufacturer, Product/Part No: as above).


2001/95/EU	General Product Safety Directive.
2014/35/EU	Low Voltage Directive (LVD).
2014/30/EU	Electromagnetic Compatibility (EMC) Directive.
2011/65/EU	Restriction on Hazardous Substances Directive (RoHS).

European harmonised Standards used for conformity assessment and declaration:

EN 61010-1: 2010 **Electrical Requirements for Laboratory Test & Measurement Equipment.**



Mark can be used only in the case of conformity assessment according to all relevant EC/EU Directives.

Dickon Knight (Director):  Date of Issue: 24th May 2023.

In the event of any correspondence concerning this product, please contact your supplying agent quoting the catalogue number and serial number shown on the apparatus rating label, together with the voltage and frequency of the mains electrical supply. This will help us to process your enquiry quickly. Any spare parts which may be required are supplied on the understanding that the replacement of those requiring the exposure of live electrical connections **will** be undertaken by an electrically qualified person.

IPC Irwin

Holker School, Cark in Cartmel, Grange over Sands, Cumbria LA11 7PQ, United Kingdom

Call: 015395 58555 **Email:** sales@ipcirwin.com