UNWIRED THINGS

IKONN TF-1

Model: FC-IR100

EN Instruction Manual

DA Brugsanvisning

NO Bruksanvisning

SV Bruksanvisning

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IS Leiðbeiningar

Manual Ver.: 1.0

Issuing Date: 03/03/2021

Introduction

Thank you for purchasing this Dual-mode Infrared Thermometer. It has been carefully developed for accurate, safe and fast temperature measurements in the ear and forehead.

Please read these instructions carefully before using this product and keep the instructions and the thermometer in a safe place.

Package Contents

No.	Name	Quantity
1	Infrared Thermometer	1
2	Pouch	1
3	Battery (AAA, optional)	2
4	User Manual	1

Contents

1.	WARNINGS AND PRECAUTIONS	4
2.	Product description	5
1)	Overview	5
2)	Structure	5
3)	Operation principle	6
4)	Indications for use	6
5)	Contraindications	6
3.	Features	6
4.	Product Structure	7
5.	Display description	7
6.	How to use your thermometer	8
1)	Take your forehead temperature	8
2)	Take your ear temperature	8
3)	Take room/object temperature	9
4)	After a measurement	9
5)	Read your Temperature	10
6)	Switching between mute and un-mute	10
7)	Checking 35 memory records	10
8)	°C/°F conversion	10
9)	Temperature compensation adjustment	10
10)	To turn off	11
11)	Replace the battery	11
7.	Temperature taking tips	11
8.	Care and cleaning	12
9.	Error and Troubleshooting	13
10.	Specifications:	14
11.	Symbols:	15
12.	EMC information	16
13.	Warranty and Customer Service	20

1. WARNINGS AND PRECAUTIONS

- 1) Keep out of reach of children under 12 years.
- 2) Never immerse the thermometer in water or other liquids (not waterproof). For cleaning and disinfection, please follow the instructions in the "Care and Cleaning" section.
- 3) Never use the thermometer for purposes other than those it has been intended for. Please follow the general safety precautions when using on children.
- 4) Keep the thermometer away from direct exposure to the sun and keep it in a dust-free, dry area, well-ventilated place at a temperature between 10 °C (50 °F) - 40 °C (104 °F). Do not use the thermometer in high humidity environments (>95% RH).
- 5) Do not use the thermometer if there are signs of damage on the measuring sensor or on the instrument itself. If damaged, do not attempt to repair the instrument! Please contact dealer.
- 6) This thermometer consists of high-quality precision parts. Do not drop the instrument. Protect it from severe impact and shock. Do not twist the instrument or the measuring sensor.
- 7) Please consult your doctor if you see symptoms such as unexplained irritability, vomiting, diarrhea, dehydration, changes in appetite or activity, seizure, muscle pain, shivering, stiff neck, pain when urinating, etc., even in the absence of fever.
- 8) Even in the absence of fever, those who exhibit a normal temperature may still need to receive medical attention. People who are on antibiotics, analgesics or antipyretics should not be assessed solely on temperature readings to determine the severity of their illness.
- 9) Temperature elevation may signal a serious illness, especially in adults who are old, frail, have a weakened immune system, or neonates and infants. Please seek professional advice immediately when there is a temperature elevation and if you are taking temperature for whom are:
 - Over 60 years of age (Fever may be blunted or even absent in elderly patients)
 - Having diabetes mellitus or a weakened immune system (e.g., HIV-positive, cancer, chemotherapy, chronic steroid treatment, splenectomy)
 - Bedridden (e.g. nursing home patient, stroke, chronic illness)
 - A transplant patient (e.g. liver, heart, lung, kidney)
- 10) This thermometer is not intended for pre-term babies or small-for-gestational age babies.

This thermometer is not intended to interpret hypothermic temperatures. Do not allow children to take their temperatures unattended.

- 11) Use of this thermometer is not intended as substitute for consultation with your physician or pediatrician. It is for household use only.
- 12) Clean the thermometer probe after each use.
- 13) Do not use the thermometer on newborns or for continuous temperature monitoring purposes.
- 14) Do not take a measurement while or immediately after nursing a baby.
- 15) Patients should not drink, eat or be physically active before/while taking the measurement.

2.Product description

Overview

Infrared Thermometer measures the body temperature based on the infrared energy emitted from the eardrum or the forehead. Users can obtain measurement results quickly after correctly positioning the temperature probe in the ear canal or on the forehead.

Normal body temperature is a range. The following tables show that this normal range also varies depending on the location. Therefore, readings from different locations should not be directly compared with each other. Tell your doctor what type of thermometer you used to measure your temperature and where on your body you measured it. Keep this in mind when self-diagnosing.


	Measurements
Forehead temperature	36.1 °C to 37.5 °C (97 °F to 99.5 °F)
Ear temperature	35.8 °C to 38 °C (96.4 °F to 100.4 °F)
Oral temperature	35.5 °C to 37.5 °C (95.9 °F to 99.5 °F)
Rectal temperature	36.6 °C to 38 °C (97.9 °F to 100.4 °F)
Axillary temperature	34.7 °C to 37.3 °C (94.5 °F to 99.1 °F)

Structure

The thermometer consists of a shell, an LCD display, a measurement button, an audible signal, an infrared temperature sensor and a microprocessor.

Operation principle

The infrared temperature sensor detects the infrared energy emitted by the eardrum or skin surface. After the energy has been focused by a lens, it is converted into a temperature value by the thermocouples and measuring circuits.

**Indications for use**

The dual-mode infrared thermometer is designed to measure human body temperature. The forehead mode is suitable for people of all ages, while the eardrum mode is suitable for people aged three months and older.


Contraindications

Do not use the thermometer if the ear is infected with otitis or suppuration.

3.Features

- Quick measurement, less than 1 second
- Accurate and reliable
- Easy to use, one-button design, for measuring the ear and forehead
- Multifunctional, can measure the temperature of the ear, forehead, room, milk, water and objects.
- 35 memory sets, easily retrievable
- Switch between mute and unmute
- Fever alarm function, indicated by orange and red lights.
- Switch between °C and °F
- Auto shut-down and power-saving


4.Product Structure

- ① Temperature light ② Measure button
 ③ LCD display screen ④ Memory/mute button
 ⑤ Probe ⑥ Probe cover (remove before measuring ear temperature)
 ⑦ Rating label ⑧ Battery cover

5.Display description

1. Temperature value
2. Object temperature mode
3. Forehead temperature mode
4. Ear temperature mode
5. Battery level
6. Mute/Un-mute icon
7. Fahrenheit / Celsius degrees
8. Ear cover
9. Memory recall


6. How to use your thermometer

When using the thermometer for the first time, please load the batteries.

1) Take your forehead temperature

Position the thermometer probe in the centre of the forehead, hold it in contact, press and hold the **Measure button**, slide it towards the temple, then release the button. A beep will sound and you can read the value.

NOTE: The forehead measurement is a guideline. The measured forehead temperature may differ from your actual body temperature by up to 1 °F/0.5 °C. Please note the factors that affect accuracy, as described in the sections "Temperature taking tips" and "WARNINGS AND PRECAUTIONS".

⚠ If the eyebrow area is covered with hair, sweat or dirt, please clean this area beforehand to improve measurement accuracy.

⚠ Always check that the lens is clean.


⚠ Always ensure that the user and the thermometer have been in the same room for at least 30 minutes before taking a measurement.

2) Take your ear temperature

Take off the probe cover, fit the probe snugly into the ear canal – Press and release the **measure button** within one second. A beep will sound and you can read the value.

⚠ Please ensure that the ear is clean and free of earwax or other foreign objects.

⚠ The reading on the right ear may differ from the reading on the left ear. Therefore, always measure the temperature in the same ear.

- ⚠ Children under 1 year of age: Pull the ear straight back.
 Children aged 1 year to adulthood: Pull the ear upwards and backwards.

⚠ Do not force the thermometer into the ear canal. Doing so may damage the ear canal.

⚠ When taking an adult's temperature, gently pull the ear upwards and backwards to ensure that the ear canal is straight so that the temperature probe can receive the infrared radiation from the eardrum.

⚠ Be careful when taking the temperature of a child whose ear canal is small.

3) Take room/object temperature

When the thermometer is turned off, press the **Memory/Mute-unmute button** for 3 seconds until it shows . Then press the **Measure button** to measure. Hold the thermometer approximately 1–5 cm away from the object. Press and release the **Measure button** within 1 second. A beep will sound and you can read the value.

4) After a measurement

After completing the measurement, remove the thermometer from the forehead/ear and read the temperature.
 After each measurement, you can enter the recall mode and view previous temperature readings.

⚠ Do not hold the thermometer in your hand for long periods of time, as it is sensitive to ambient temperature.

⚠ Clean the thermometer with a soft cloth after each measurement and store it in a dry, well-ventilated place.

⚠ You should wait at least 10 seconds between measurements.

⚠ It is dangerous to self-diagnose or self-treat based on the measurement results obtained. Please consult a doctor for this.

5) Read your Temperature

T displays a temperature value. In forehead or ear mode.


1. When $32\text{ }^{\circ}\text{C}$ ($89.6\text{ }^{\circ}\text{F}$) $\leq T \leq 37.3\text{ }^{\circ}\text{C}$ ($89.6\text{ }^{\circ}\text{F}$) $\leq T \leq 99.2\text{ }^{\circ}\text{F}$, the green LED lights up for 3 seconds and a long beep sounds.

2. If $37.4\text{ }^{\circ}\text{C} \leq T \leq 37.9\text{ }^{\circ}\text{C}$ ($99.3\text{ }^{\circ}\text{F} \leq T \leq 100.3\text{ }^{\circ}\text{F}$), the orange lamp will light up for 3 seconds, 3 short beeps will sound and the value on the LCD display will flash. This is a warning that you may have a slight fever.

3. If $38\text{ }^{\circ}\text{C} \leq T \leq 42.9\text{ }^{\circ}\text{C}$ ($100.4\text{ }^{\circ}\text{F} \leq T \leq 109.2\text{ }^{\circ}\text{F}$), the red lamp will light up for 3 seconds, 5 short beeps will sound, and the value on the LCD display will flicker. This is a warning that you may have a high fever.

6) Switching between mute and un-mute

When the thermometer is turned on, press and hold the **Memory/Mute-unmute button** for 2-3 seconds to switch from mute to unmute.

Checking 35 memory records

When the thermometer is turned on or off, press the **Memory/Mute** button to enter memory mode, and press this button again to view the 35 memory records in sequence. If no value is recorded, "—M" will be displayed.


°C/°F conversion

Open the battery cover, use the toggle switch to change the °C/°F.

Temperature compensation adjustment

When the thermometer is switched on, press both the **Memory/Mute-unmute button** and the **Measure button** simultaneously for 2–3 seconds to enter temperature compensation mode. By pressing the **Memory/Mute-unmute button**, you can adjust the temperature from ± 0.0 to ± 2.0 .

Note: All the future temperature you are taking will be automatically added the value you are adjusted.

To turn off

The device switches off automatically after 10 seconds of inactivity. Alternatively, you can press and hold the **Measure button** for 5–7 seconds.

Caution

1. All memory records will be lost if the battery is removed or reinserted.
2. All settings will be reset to default when the battery is removed. If you need to adjust the settings, switch on the device and make the new settings.

Replace the battery.

Slide the battery cover off in the direction indicated. Insert two AAA batteries correctly into the compartment.

 Remove the batteries if the thermometer will not be used for more than two months.

7. Temperature taking tips

1) It is important to know each individual's normal body temperature when they are healthy. This is the only way to accurately diagnose fever. Take measurements twice a day (early in the morning and late in the afternoon). Take the average of the two temperatures to calculate the normal oral equivalent temperature. Always measure the temperature in the same place, as temperature readings can vary at different points on the forehead.

2) A child's normal body temperature can be up to 37.7 °C (99.9 °F) or up to 36.11 °C (97.0 °F). Please note that this device displays 0.5 °C (0.9 °F) less than a digital rectal thermometer.

- 3) External factors can influence ear temperature, including:
- lying on one ear or the other
 - covering your ears
 - being exposed to very hot or very cold temperatures, or
 - recently swimming or bathing


In these cases, remove the person from the situation and wait 20 minutes before taking their temperature.

Use the untreated ear if prescription ear drops or other ear medications have been placed in the ear canal.


- 4) Holding the thermometer in your hand for too long before taking the measurement may cause the device to warm up. This means that the measurement may be inaccurate.
- 5) Patients and thermometers should remain in a room with stable conditions for at least 30 minutes.
- 6) Before placing the thermometer sensor on the forehead, remove any dirt, hair or sweat from the forehead area. Wait 10 minutes after cleaning before taking the measurement.
- 7) Carefully clean the sensor with an alcohol swab and wait 5 minutes before taking a measurement on another patient. Wiping the forehead with a warm or cold cloth may affect your readings. It is advised to wait 10 minutes before taking a reading.
- 8) In the following situations, it is recommended to take 3-5 temperatures at the same spot and use the highest value as the reading:
- Newborns in the first 100 days.
 - Children under three years of age with a weakened immune system, where the presence or absence of fever is critical.
 - When the user is using the thermometer for the first time, until they have familiarised themselves with the device and are obtaining consistent readings.

8. Care and cleaning

Clean the thermometer housing and measuring probe with an alcohol swab or a cotton swab moistened with 70% alcohol. Once the alcohol has dried completely, you can take a new measurement. Ensure that no liquid gets inside the thermometer. Never use abrasive cleaning agents, thinners or benzene for cleaning, and never immerse the device in water or other cleaning fluids. Be careful not to scratch the surface of the LCD screen.

9.Error and Troubleshooting

Symptom	Possible cause	Description & Solution
The device cannot be switched on.	The battery level is too low.	Replace with a new battery.
	The polarities of the batteries are reversed.	Ensure that the batteries are inserted correctly.
	The thermometer is damaged.	Contact your dealer
The measured value is too low	The probe lens is dirty.	Clean the lens with a cotton swab.
	The distance between the object and the target is too great	Keep the thermometer in contact with the forehead or insert the probe into the ear canal.
	You have just come from a cold environment	Stay in a warmer room for at least 30 minutes in a warmer room before taking a measurement
The reading is too high	You have just come from a hot environment.	Stay in a sufficiently cool room for at least 30 minutes before taking a measurement.
	The ambient temperature is not within the permissible range.	3 short beeps and red backlight for 3 seconds. Take the measurement at an ambient temperature between 10 °C (50.0 °F) and 40 °C (104 °F).
	Memory error	3 short beeps and red backlight for 3 seconds. Contact your retailer.
	In ear or forehead mode, $T > 42.5\text{ °C}$ (108.6 °F)	3 short beeps and red backlight for 3 seconds.

Symptom	Possible cause	Description & solution
	In ear or forehead mode, $T < 32\text{ }^{\circ}\text{C}$ ($89.6\text{ }^{\circ}\text{F}$)	3 short beeps and 3 seconds of red backlight.
	$2.5\text{ V} \pm 3\% \leq \text{power voltage} \leq 2.6\text{ V} \pm 3\%$	The battery level is low. It is recommended that you replace the battery, but you can continue to use the device.
	The supply voltage is below $2.5\text{ V} \pm 3\%$.	It will switch off automatically after 30 seconds. Please replace the battery with a new one.

10.Specifications:

Product name	Dual-mode infrared thermometer	
Power	DC1.5 V x 2	
Measuring range	Ear and forehead: $32.0\text{ }^{\circ}\text{C} - 42.9\text{ }^{\circ}\text{C}$ ($89.6\text{ }^{\circ}\text{F} - 109.2\text{ }^{\circ}\text{F}$)	
	Object: $0\text{ }^{\circ}\text{C} - 10\text{ }^{\circ}\text{C}$ ($32\text{ }^{\circ}\text{F} - 212\text{ }^{\circ}\text{F}$)	
Accuracy (laboratory)	Ear and forehead mode	$\pm 0.2\text{ }^{\circ}\text{C}/\pm 0.4\text{ }^{\circ}\text{F}$
	Object mode	$\pm 1.0\text{ }^{\circ}\text{C}/1.8\text{ }^{\circ}\text{F}$
Screen resolution	$0.1\text{ }^{\circ}\text{C}/^{\circ}\text{F}$	
Automatic switch-off	$10\text{s} \pm 1\text{s}$	
Memory	35 temperature measurement groups.	
Operating conditions	Temperature: $10^{\circ}\text{C}-40^{\circ}\text{C}$ ($50^{\circ}\text{F}-104^{\circ}\text{F}$)/ Humidity: 15–95% relative humidity, non-condensing Air pressure: 86-106 kPa	
Battery	2*AAA, can be used more than 3000 times	
Weight & Dimensions	66 g (without battery), 163.3 x 39.2 x 38.9 mm	

11.Symbols:

Symbol	Description
	Type BF applied part
	MANUFACTURER
	Refer to instruction manual/booklet
	The marking of electrical and electronics devices according to Directive 2012/19/EC. The device accessories and the packaging have to be disposed of correctly at the end of the usage. Please follow Local Ordinances or Regulations for disposal
	SERIAL NUMBER
	Batch number
	IMPORTANT Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.
IP22	2 Protected against solid foreign objects of 12,5 mm Ø and greater ; 2 If keep the thermometer is at a 15 degree angle, it still protected against dripping water.
	AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
	Transport package shall be kept away from rain.
	Keep away from sunlight.
	Fragile, handle with care
	Unique device identifier
	Date of manufacture
	Model number
	Importer
	Made in China
	Medical device
	CE mark

12. EMC information

IEC 60601-1-2:2014/AMD1:2020 ME EQUIPMENT and ME SYSTEMS identification, marking and documents for Class B product

Instructions for use

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments and so on.

Warning: Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

If any: A list of all cables and maximum lengths of cables (if applicable), transducers and other ACCESSORIES that are replaceable by the RESPONSIBLE ORGANIZATION and that are likely to affect compliance of the ME EQUIPMENT or ME SYSTEM with the requirements of Clause 7 (EMISSIONS) and Clause 8 (IMMUNITY). ACCESSORIES may be specified either generically (e.g. shielded cable, load impedance) or specifically (e.g. by MANUFACTURER and EQUIPMENT OR TYPE REFERENCE).

If any: The performance of the ME EQUIPMENT or ME SYSTEM that was determined to be ESSENTIAL PERFORMANCE and a description of what the OPERATOR can expect if the ESSENTIAL PERFORMANCE is lost or degraded due to EM DISTURBANCES (the defined term "ESSENTIAL PERFORMANCE" need not be used).

Technical description

1. all necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.

2. Guidance and manufacturer's declaration -electromagnetic emissions and Immunity.

Table 1

Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	N/A
Voltage fluctuations/ flicker emissions IEC 61000-3-3	N/A

Table 2

Guidance and manufacturer's declaration - electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/ burst IEC 61000-4-4	Power supply lines: ± 2 kV 100 kHz repetition frequency	N/A
Surge IEC 61000-4-5	line(s) to line(s): ± 1 kV.	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 25°, 270° and 315° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0% 300 cycle	N/A
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz

Conducted RF * IEC61000-4-6	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	N/A
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
Proximity magnetic fields IEC 61000-4-39	30 kHz: 8A/m 134.2 kHz: 65A/m 13.56 MHz: 7.5A/m	30 kHz: 8A/m 134.2 kHz: 65A/m 13.56 MHz: 7.5A/m
NOTE UT is the a.c. mains voltage prior to application of the test level. Remark:“*” : Only test for Patient Coupling line.		

Table 3

Guidance and manufacturer's declaration - electromagnetic Immunity		
Test frequency	Modulation	IMMUNITY TEST LEVEL (A/m)
30 kHz	CW	8
134,2 kHz	Pulse modulation ^a 2,1 kHz	65 ^b
13,56 MHz	Pulse modulation ^a 50 kHz	7,5 ^b
a) The carrier shall be modulated using a 50% duty cycle square wave signal. b) r.m.s., before modulation is applied.		

Table 4

Guidance and manufacturer's declaration - electromagnetic Immunity								
	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	IEC 60601-1-2 Test level (V/m)	Compliance level (V/m)
Radiated RF IEC610-00-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380–390	TETRA 400	Pulse modulation 18 Hz	1,8	0.3	27	27
	450	430–470	GMS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28	28
	710	704–787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0.3	9	9
	745							
	780							
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDM A 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28	28
	870							
	930							
	1720	1700 – 1990	GSM 1800; CDM A 1900; GSM 1900; DECT ;	Pulse modulation 217 Hz	2	0.3	28	28
	1845							

	1970		LTE Band 1, 3, 4, 25; UMTS					
	2450	2400 – 2570	Bluetooth, WLAN , 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28	28
	5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9	9
	5500							
	5785							

13. Warranty and Customer Service

The device is under warranty for 24 months from the date of purchase. The batteries, the packaging and any damage caused by improper use are not covered by the warranty.

Excluding the following user-caused failures:

- 1.Failure resulting from unauthorized disassembly and modification.
- 2.Failure resulting from an unexpected dropping during application or transportation.
- 3.Failure resulting from not following the instructions in the operating manual.

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CTI U.S. Inc.
Suite 230,1455 Cti U.S. Inc.
Shenzhen Finicare Co., Ltd